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
# NKV-550 Series Ventilator System Specifications

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### **Manufacturer Address**

Nihon Kohden OrangeMed, Inc.  
1800 E. Wilshire Avenue  
Santa Ana, CA 92705  
United States of America

### **Website**

[www.orange-med.com](http://www.orange-med.com)  
[www.nihonkohden.com](http://www.nihonkohden.com)

### **Email**

Customer Service Department: [customers@orange-med.com](mailto:customers@orange-med.com)  
Technical Service and Support: [techservice@orange-med.com](mailto:techservice@orange-med.com)  
Clinical Education and Support: [clinical@orange-med.com](mailto:clinical@orange-med.com)

### **Telephone**

(Monday through Friday, 8:00 am to 5:00 pm, Pacific Time)

### **Main Office: +1 (949) 502-6448**

Customer Service: Ext: 7101  
Technical Support: Ext. 7102  
Clinical Support: Ext. 7103

### **EC Representative**

NIHON KOHDEN EUROPE GmbH  
Raiffeisenstrasse 10  
D-61191 Rosbach  
Germany  
Tel +49 (6003) 827-0  
Fax +49 (6003) 827-599

### **Copyright Information**

© Nihon Kohden OrangeMed, Inc. All rights reserved. Nihon Kohden OrangeMed, Inc.'s NKV-550 Ventilator System is manufactured in accordance with Nihon Kohden OrangeMed, Inc.'s proprietary information.

Patent pending.

The information in this manual is the sole property of Nihon Kohden OrangeMed, Inc. and may not be duplicated without permission. This manual may be revised or replaced by Nihon Kohden OrangeMed, Inc. at any time and without notice.

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# Specifications

## Configurations

The NKV-550 Series Ventilator System configurations are detailed in Table 1-1.

Table 1-1: NKV-550 Series Feature Comparison Chart

Feature	NKV-550-N (Neonatal)	NKV-550-S (Standard)	NKV-550-U (Universal)
Neonatal Patient Type	●		●
Pediatric Patient Type 3.1	●	●	●
Adult Patient Type		●	●
Extended Battery	■	■	●
Neonatal On-airway Flow Sensor	●		●
Aerogen Nebulizer	■	■	■
Nihon Kohden SpO <sub>2</sub> *	●	●	●
Nihon Kohden CO <sub>2</sub> *	●	●	●

\* Sensor purchased separately

● = Standard Feature

■ = Optional Feature

## Patient Type

- Adult
- Pediatric
- Neonate

## Patient Data

ID	0 - 9999
Gender	Male / Female
Height	130 cm - 200 cm
Body Weight (BW)	0.30 kg or higher
Predicted Body Weight (PBW)	Calculated from Gender and Height Inputs

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**Ventilation Modes**

4.1

Invasive Ventilation	3.1 4.4	A/CMV-PC SIMV-PC-PS SPONT-CPAP APRV	A/CMV-VC SIMV-VC-PS SPONT-PS	A/CMV-PRVC SIMV-PRVC-PS SPONT-VS
Non-invasive Ventilation	4.4	A/CMV-PC SIMV-PC-PS SPONT-CPAP APRV nCPAP	SPONT-PS	
Oxygen Therapy	4.4	O <sub>2</sub> Therapy		

**Ventilation Settings**

Tidal Volume (VT)	2 to 100 mL 5 to 100 mL 20 to 1000 mL 100 to 3000 mL	Neonate (PRVC/VS) Neonate Pediatric Adult	4.8
Pressure Control (P <sub>INSP</sub> or ΔPC)	2 to 60 cmH <sub>2</sub> O 2 to 70 cmH <sub>2</sub> O 2 to 80 cmH <sub>2</sub> O	Neonate (60-PEEP) Pediatric (70-PEEP) Adult (80-PEEP)	4.12
Pressure Support (PS)	0 to 60 cmH <sub>2</sub> O 0 to 70 cmH <sub>2</sub> O 0 to 80 cmH <sub>2</sub> O	Neonate (60-PEEP) Pediatric (70-PEEP) Adult (80-PEEP)	4.13
PEEP	0 to 30 cmH <sub>2</sub> O 0 to 40 cmH <sub>2</sub> O 0 to 50 cmH <sub>2</sub> O	Neonate Pediatric Adult	4.14
CPAP	0 to 30 cmH <sub>2</sub> O 0 to 40 cmH <sub>2</sub> O 0 to 50 cmH <sub>2</sub> O	Neonate Pediatric Adult	
P <sub>HIGH</sub>	1 to 50 cmH <sub>2</sub> O		
P <sub>LOW</sub>	0 to 49 cmH <sub>2</sub> O		
T <sub>HIGH</sub>	0.1 to 30 sec		
T <sub>LOW</sub>	0.1 to 30 sec		

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NKV-550 Series Ventilator System | 5

## Specifications

Flow Type	Square, Descending 50%	4.5
4.10 Flow Rate (Flow)	<p>Volume Control:</p> <p>1 to 30 L/min      Neonate          1 to 60 L/min      Pediatric          1 to 150 L/min    Adult</p> <p>PC, PS, PRVC, VS, Spont:          Up to 180 L/min</p>	
	<p>O<sub>2</sub> Therapy:</p> <p>OFF, 1 to 15 L/min      Neonate          OFF, 1 to 30 L/min      Pediatric          OFF, 1 to 60 L/min    Adult</p>	
Inspiratory Pause (Pause)	OFF, 0.1 to 2.0 sec	4.16
Inspiratory Time (Ti)	<p>0.20 to 3.0 sec      Neonate / Pediatric (VC)          0.20 to 5.0 sec      Adult (VC)          0.20 to 10.0 sec    All patient sizes (PC)</p>	4.15
I:E Ratio (I:E)	4.0:1 to 1:299	
Respiratory Rate (RR)	<p>1 to 150 bpm      Neonate          1 to 120 bpm      Pediatric          1 to 80 bpm      Adult</p>	4.9
Oxygen % (FiO <sub>2</sub> )	21 to 100%	
Trigger Type (P <sub>TRIG</sub> or F <sub>TRIG</sub> ):	<p>Pressure Trigger:          0.1 to 20 cmH<sub>2</sub>O</p> <p>Flow Trigger:          0.1 to 20 L/min    Adult          0.1 to 15 L/min    Pediatric          0.1 to 10 L/min    Neonate</p>	4.11
Slope	5% (Slowest) to 100% (Fastest)	4.17
Expiratory Trigger (ET%)	1 to 80%	4.18
Maximum Inspiratory Time of PS (TiMax PS)	<p>0.3 to 1.0 sec      Neonate          0.5 to 1.5 sec      Pediatric          0.8 to 2.0 sec      Adult</p>	
Sigh	<p>OFF, ON          Factor:      1.1 to 1.5          Interval:    30 to 100 (control breaths)</p>	
Apnea Ventilation	<p>OFF, ON          VT apn          RR apn</p>	

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Tube Compensation (Tube Comp)	4.6	OFF, ON Tube type: ETT or Trach Tube ID: 2.0 to 10.0 mm Comp%: 0 to 100%
Leak Compensation	4.7	Invasive: ON/OFF up to 10 L/min      Neonate up to 15 L/min      Pediatric up to 25 L/min      Adult  Non-invasive: ON only up to 15 L/min      Neonate up to 40 L/min      Pediatric up to 65 L/min      Adult  Max Vol LC (VC only) 0 - 50 mL      Neonate 0 - 100% of set VT      Pediatric / Adult

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## Specifications

Driving Pressure ( $P_{DRIVING}$ )	0 - 99.9 cmH <sub>2</sub> O
Driving Pressure, Estimated ( $P_{DRIVING-EST}$ )	0 - 99.9 cmH <sub>2</sub> O
5.2.1    Inspiratory Tidal Volume (VT <sub>i</sub> )	0 - 3,500 mL
Tidal Volume (VT)	0 - 3,500 mL
Tidal Volume per Kg (VT/kg)	0 - 50 mL/kg
5.2.2    Minute Volume (MV)	0.00 - 99.9 L
Spontaneous Minute Volume (MV <sub>SPONT</sub> )	0.00 - 99.9 L
Leak at PEEP	0 to 200 L/min
Leak Volume % (Leak %)	0 - 100%
5.2.3    Leak Volume (V <sub>LEAK</sub> )	0 - 3,000 mL
Total Respiratory Rate (RR <sub>TOT</sub> )	0 - 200 bpm
5.3.1    Spontaneous Respiratory Rate (RR <sub>SPONT</sub> )	0 - 150 bpm
5.3.2    Mandatory I:E Ratio (I:E)	16.0:1 to 1:299
APRV T <sub>H</sub> and T <sub>L</sub> Ratio (T <sub>H</sub> :T <sub>L</sub> )	150:1 to 1:150
5.3.2    Spontaneous Inspiratory Time (T <sub>ISPONT</sub> )	0.10 - 9.99 s
5.3.2    Spontaneous Duty Cycle (T <sub>i</sub> /T <sub>TOT</sub> )	10 - 90%
5.4.1    Static Inspiratory Resistance (R <sub>i-STAT</sub> )	1 - 200 cmH <sub>2</sub> O/L/s
5.4.2    Static Compliance (C <sub>STAT</sub> )	0.1 - 120 mL/cmH <sub>2</sub> O
Static Compliance per kg (C <sub>STAT/kg</sub> )	0.00 - 5 mL/cmH <sub>2</sub> O/kg
Expiratory Resistance (R <sub>E</sub> )	1 - 200 cmH <sub>2</sub> O/L/s
Dynamic Resistance, Estimated (R <sub>EST</sub> )	1 - 200 cmH <sub>2</sub> O/L/s
5.4.2    Dynamic Compliance, Estimated (C <sub>EST</sub> )	0.1 - 120 mL/cmH <sub>2</sub> O
5.3.2    Exhalation Time Constant (TC <sub>E</sub> )	.01 - 6 s
Imposed Work of Breathing (WOB <sub>IMP</sub> )	0 - 99.9 J/min
C20/C	0.1 - 3.0

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## Specifications

5.4.3	Rapid Shallow Breathing Index (RSBI)	0 - 9999 bpm/L
	Rapid Shallow Breathing Index per kg (RSBI/kg)	0 - 300 bpm/mL/kg
	Oxygen Concentration FiO <sub>2</sub>	18 – 100%
	Oxygen Pulse Saturation (SpO <sub>2</sub> ) and SQI bar graph	0 - 100%
	Pulse Rate (PR)	30 – 300 bpm
	Pulse-Amplitude Index (PI)	0.01 - 100 %
	End tidal CO <sub>2</sub> (EtCO <sub>2</sub> )	0 to 150 mmHg
Waveforms and Loops		Pressure Waveform Flow Waveform Volume Waveform  Auxiliary Pressure Waveform  Pressure-Volume Loop Flow-Volume Loop
		3.3

### Alarm Audio Pause

Audio Pause	2 min (max)
-------------	-------------

### Adjustable Alarms

Airway Pressure (Paw), High	5 - 100 cmH <sub>2</sub> O 5 - 120 cmH <sub>2</sub> O	Neo/Ped Adult	5.6.1
Minute Ventilation (MV), High	0.02 - 20.0 L	Neonate	5.6.2
	0.03 - 40.0 L 0.03 - 60.0 L	Pediatric Adult	
Minute Ventilation (MV), Low	0.01 - 19.0 L 0.02 - 39.0 L 0.02 - 59.0 L	Neonate Pediatric Adult	
	OFF available, NIV only		
Tidal Volume mL/kg (VT), High	2 - 30 L, OFF		5.6.3
Tidal Volume mL/kg (VT), Low	OFF, 1 - 29 L		

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## Specifications

Respiratory Rate (RR), High	10 - 150 bpm, OFF 10 - 120 bpm, OFF	Neo/Ped Adult	5.6.4
Apnea	5 to 60 sec  OFF available, NIV only		5.7
Leak, High	20 - 95%, OFF		
SpO <sub>2</sub> , High	51 - 100%, OFF		
SpO <sub>2</sub> , Low	OFF, 50 - 99%		
Pulse Rate (PR), High	31 - 300 bpm, OFF		
Pulse Rate (PR), Low	OFF, 30 - 299 bpm		
EtCO <sub>2</sub> , High	2 - 99 mmHg, OFF		
EtCO <sub>2</sub> , Low	OFF, 1 - 98 mmHg		

## Quick Access Buttons

Home	Inspiratory Hold
Panel Lock	Expiratory Hold
Elevated O <sub>2</sub>	Screen Brightness
Manual Breath	Help

## Applications (Apps)

Standby	Data Retrieval
Open Airway Suctioning	Camera
In-line Airway Suctioning	Sensors
NIF/PiMax Maneuver	Custom Settings
P0.1 Measurement	Video
Low Flow PV Maneuver	
Volumetric Capnography	<b>Optional Apps</b>
Spontaneous Breathing Trial	Recruitability Assessment
Auxiliary Pressure	Recruitment Maneuver
Trends	5.5 PEEP Titration
Logs	Transpulmonary Pressure

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## Institutional Settings

System Settings	Language Pressure units Patient height units Communication protocol Nurse call settings
-----------------	---

## Gases

O <sub>2</sub> Supply	Input: 25 to 87 psi (172 to 600 kPa) Vmax: 180 L/min
Air Supply	Input: 25 to 87 psi (172 to 600 kPa) Vmax: 180 L/min

## Physical Specifications

Dimensions	Display (not including mount)
	Height            46.7 cm (18 25/64 in) Width            34.4 cm (13 1/2 in) Depth            58.0 cm (2 9/32 in)
4.3	Breath Delivery Unit
	Height            27.3 cm (10 3/4 in) Width            43.5 cm (17 1/8 in) Depth            48.4 cm (19 in)
Weight	Standard Configuration (Display mounted on BDU)
	Height            69.8 cm (27 1/2 in) Max. tilt Width            43.5 cm (17 1/8 in) Depth            53.5 cm (21 in) Max. tilt
Weight	Display (not including mount)
	Weight            5.2 kg (11.4 lbs)
Weight	Breath Delivery Unit
	Weight            18 kg (39.7 lbs)
Weight	Standard Configuration (Display mounted on BDU)
	Weight            24.2 kg (53.4 lbs)

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## Specifications

## Environmental

Operation	Temperature: 10°C to 40°C (50°F to 104°F) Humidity: 10% to 95% non-condensing Atmospheric pressure: 700 hPa to 1060 hPa (10.15 psi to 15.37 psi) Altitude: -411.5 m to 3048 m (-1350 ft to 10000 ft)
Storage	Temperature: -20°C to 50°C (-68°F to 122°F) Humidity: 10 to 95% non-condensing Atmospheric pressure: 500 hPa to 1060 hPa (7.25 psi to 15.37 psi) Altitude: 6096 m max (20000 ft max)

## Powers

Mains (AC Power)	3.4	100 to 240 V, 50/60 Hz Power consumption: 50 - 200 W Fuse 1 & 2: T3.15A, (250V), H
Extended Battery		Li-ion, 14.4 V, 9.0 Ah Run time: 2 hr. 10 min (new and fully charged)
Backup Battery	3.4	Li-ion, 14.4 V, 4.1 Ah Run time: 50 min. (new and fully charged)

The NKV-550 ventilator is classified as IP21. This IEC ingress protection classification indicates protection by the enclosure against ingress of foreign objects and dripping liquid.

The NKV-550 ventilator has been tested to be compliant with the following standards:

Standard	Testing
ISO 80601-2-12 2011/04/15 Ed: 1; CORR 1: 2011/10/15	Medical Electrical Equipment - Part 2-12: Particular Requirements For Basic Safety And Essential Performance Of Critical Care Ventilators
ANSI AAMI ES60601-1 2005+AC1; A2 (R2012)	Medical Electrical Equipment - Part 1: General Requirements For Basic Safety And Essential Performance
IEC 60601-1-2 2014/02/25 Ed: 4	3.6 Medical Electrical Equipment - Part 1-2: General Requirements for Safety - Collateral Standard: Electromagnetic Compatibility - Requirements and Tests
IEC 60601-1-6 2013/10/28 Ed: 3:1	Medical Electrical Equipment - Part 1-6: General requirements for safety - Collateral Standard: Usability
IEC 60601-1-8 2012/11/28 Ed: 2.1	Medical Elec. Equip. - Part 1-8: General Req. for Basic Safety & Essential Perf. - Collateral Standard: General Req., Tests & Guidance for Alarm Systems in Medical Elec. Equip. & Medical Elec. Systems

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Specifications

Standard	Testing
IEC 62133-2 2017/02/07 Ed: 1.0	Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary lithium cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems
ISO 80601-2-55 2011/12/15 Ed:1	Medical Electrical Equipment - Part 2-55: Particular Requirements For The Basic Safety And Essential Performance Of Respiratory Gas Monitor
ISO 10993-1 2009/10/15 Ed: 4; TC 1 2010	Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing within a Risk Management Process
ISO 18562-1 to 3 2017/03 Ed: 1	Biocompatibility Evaluation of Breathing Gas Pathways In Healthcare Applications

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
## Chapter 3: System Overview

### System Overview

The following chapter provides an introduction to the NKV-550 ventilator system, including the breath delivery unit, graphic user interface, touchscreen, key components, connectors, symbols, icons, and ventilator mounting options.

4.3

Figure 3-1: Main System Components Front and Back

#### 1 Graphic user interface 3.2

Allows clinicians to:

- Set ventilator control parameters, such as tidal volume and inspiratory pressure
- Set alarm limits, such as high inspiratory pressure alarm
- View monitored numeric values
- View waveforms and loops
- Operate various features via Apps

#### 2 Breath delivery unit

- Receives inputs from the electronic system and controls gas delivery to the patient
- Provides various alarms, a safety valve, and other design features to maximize patient safety

## Second Graphic User Interface 3.2

The NKV-550 ventilator is compatible with an optional NKV-550 second graphic user interface (GUI) that can be used in addition to the primary GUI. This second GUI may provide convenience in critical care environments where the primary GUI may be in a location that is difficult for clinicians to access.

The second GUI and connecting cable must be purchased additionally.


**WARNING:**

The patient has to be within the sight of clinicians who operate the second graphic user interface (GUI).

**To connect a second graphic user interface:**

- 1 Connect one end of the connecting cable to the second GUI port that is located on the back of the breath delivery unit.
- 2 Connect the opposite end of the connecting cable to the connector on the second GUI.
- 3 Use the illustration below as a guide to connecting the second GUI.

*Figure 7-8: Connection of the Second Graphic User Interface (GUI)*

**Note:**

Always ensure that the primary graphic user interface (GUI) cable is connected prior to powering on the NKV-550 ventilator.

## Chapter 5: Theory of Operations

Regulator. If circuit pressure unexpectedly increases due to an obstruction, any pressure exceeding 125 cmH<sub>2</sub>O will be "vented" through the Safety Valve relief port.

If a Device Alert occurs and the ventilator stops delivering breaths, or if the expiratory limb of the breathing circuit is occluded, or if the exhalation valve is stuck closed, the Safety Valve Solenoid activates, interrupting the back pressure holding the Safety Valve Diaphragm closed. With the back pressure now removed, the diaphragm is pushed open by the Safety Valve Spring, allowing circuit pressure to drop to ambient pressure by "venting" through the Safety Valve Relief port.

### **Expiratory Module 4.2**

When the patient exhales, the exhaled gas travels through the expiratory limb of the breathing circuit and enters the expiratory port of the ventilator where it passes through the exhalation valve and hot-wire expiratory flow sensor, then finally exits through the ventilator's exhaust port.

The Expiratory Module is heated to minimize any "rain out" from the cooling of the exhaled gas. The hot-wire flow sensor continuously measures the exhaled flow rate for the ventilator to compute the volume and flow based parameter data.

### **Exhalation Control Module**

The exhalation control module provides the precise control of the exhalation valve for both inspiration and expiration.

Driving pressure to the exhalation valve is provided by the high pressure air supply, is controlled by the Exhalation Regulator and Exhalation Proportional Valve, and monitored by the P3 transducer. Exhalation circuit pressure is monitored by the P2 transducer. An in-line filter helps protect the sensor from contamination and a purge solenoid helps to ensure that this pressure tube to the transducer remains free of moisture.

If the high pressure air supply fails to the ventilator, a crossover solenoid will activate, and the high pressure oxygen will supply the pressure to the expiratory module.

### **Auxiliary Pressure Module**

The Auxiliary Pressure Module provides an external connector to monitor an additional pressure. The P7 transducer monitors this pressure and an in-line filter helps protect the sensor from contamination.

## Chapter 6: Applications (Apps)

### After completion of the **NIF/PiMax** Maneuver:

- Normal ventilation will resume.
- A vertical cursor will be displayed on the Paw waveform. A number showing the value of the pressure at the cursor will be displayed.
- Clinicians can adjust the cursor by the encoder knob to the lowest airway pressure to obtain the **NIF/PiMax** measurement.
- Press **Accept** button to save the measured value or **Reject** button to reject the measured value.

### Spontaneous Breathing Trial (SBT)

SBT app provides a group of parameter trend graphs and real-time monitored parameters for clinicians to evaluate a patient's spontaneous breathing.

#### Applicable modes/breath types

SIMV and SPONT invasive modes and breath types.

#### Range:

Up to 120 min. (max)

#### How to Use:

- Select the **SBT** app button.
- The app dialog box will open displaying:
  - 4 parameter/trend graphs
  - 6 real-time monitored parameters
- Data points will plot each minute
- A cursor will be positioned at the left side of each graph.
  - Rotating the encoder knob will move the cursor between data points for the 4 graphs.
  - The numeric value for the cursor's location will be displayed at the top of each graph.
- Pressing a trend graph parameter cell allows for selection of alternative trend parameters.
- A *download* button allows downloading of the SBT result via a USB.

#### To perform a SBT:

- Press **START RECORDING** button to begin SBT.
- Start time of SBT and the SBT total time will be displayed.
- Press **STOP RECORDING** button to end the SBT.

### Camera

5.8

Camera app can be used to take a screenshot of the ventilator's screen

#### Applicable modes/breath types

Any time when the ventilator is ON.

#### How To Use:

- Select the Camera app button.

## Chapter 6: Applications (Apps)

- A time bar will display both the total time and elapsed time of the *PEEP-T*. Numeric values will also be displayed for elapsed time and remaining time of the *PEEP-T*.
- A PEEP time-based graph will be displayed throughout the *PEEP-T*.
- A second time-based graph will also be displayed throughout the *PEEP-T* that is clinician selectable. Default parameter is estimated dynamic compliance ( $C_{EST}$ ).
- Pressing the *END* button will cancel the *PEEP-T* and no results will be displayed.

Note:

During the *PEEP-T*, patient triggering will be disabled. If a patient effort is detected, a *Patient Effort Detected* alarm message will be displayed.

### After completion of the PEEP Titration:

- When the *PEEP-T* is completed, the ventilator will resume normal ventilation with the current PEEP control setting.
- The following results of the *PEEP-T* will be displayed:
  - *PEEP at highest compliance*: the measured PEEP at the highest compliance level.
- A vertical cursor will be displayed on both graphs along with the measured values for the current cursor position. Clinicians can move the cursor to view measured results throughout the *PEEP-T*.
- Clinician can choose to *END* the *PEEP-T* or *Go to Recruitment Maneuver*.
  - To end *PEEP-T*, press *END* button.
  - To proceed to *RM* app, press *Go to Recruitment Maneuver* button.

## Video

Video app can be used to take video of the ventilator's screen.

**Applicable modes/breath types**

Any time when the ventilator is ON.

### How To Use:

- Select the *Video* app button.
- A 30 second video clip of the ventilator screen will be recorded.
- Saved video can be retrieved through the *Data Retrieval* app. See *Data Retrieval* app for more information.

## Chapter 6: Applications (Apps)

the clinician selects "new patient", then *Custom Settings* app. The ventilator must remain in *Start Up* screen to load a new *Custom Setting Profile*.

- **Delete:** Allows for deletion of a saved *Custom Setting Profile*. Clinician selects the profile for deletion and then selects the **Delete** button.

## Trends 5.9

*Trends* app can be used to review patient ventilation parameter trend graphs for up to 168 hrs.

#### **Applicable modes/breath types**

All therapy types, modes and breath types.

## Range

**Range**  
1 to 168 hours

#### **How To Use:**

- Select the **Trends** app button.
  - Press the *Channel Selection* button to select between 1 to 4 channels for display.
  - Press the *Group* button to select one of four customizable trend parameter groupings.
  - To customize a *Group*:
    - Select the trend *Group* to be customized.
    - Select the parameter cell within the *Group* to be changed.
    - Select the replacement parameter.
  - Press the + or - buttons to adjust the trend time scale.
  - A *download* button allows downloading of the *Trends* via a USB port.

## Logs

Logs app can be used to review logs of various types of events such as ventilator setting changes and alarm violations.

#### **Applicable modes/breath types**

All therapy types, modes and breath types.

## How To Use:

- Select the **Logs** app button.
  - The app dialogue box will open and the most recent logged events will be displayed.
    - Arrow buttons allow for scrolling through the logged events.
    - *Date/Time* buttons allow for selecting a specific timeframe for review.
  - Primary ventilator control settings and User alarm settings will be listed for the highlighted event selected.
  - Logs can be filtered by either *event* or *type*.
    - To filter by *event*, select a logged event, then press the *Event* button.
    - To filter by *type*, select the logged type, then press the *Type* button.

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## Chapter 6: Applications (Apps)

- button to toggle button to *Enabled*.
- Press the **Ptp** app button.
  - The app dialogue box will open and will display the PTP waveform and the measured values of the transpulmonary pressures.
  - Pressure (Ptp) is calculated by (Paw - Paux). Measured values include *Ptp* at peak, plateau, mean and end-exp and  $\Delta$  pressure.
  - Zeroing of the *Auxiliary* pressure transducer can be performed in the *Sensors* app.

## Low Flow P/V Maneuver 5.5

*Low flow P/V maneuver* app can be used to perform a pressure-volume (P/V) loop with a constant flow at a low flow rate. It allows clinicians to measure respiratory compliance, lower inflection point (LIP) and upper inflection point (UIP) from the low flow pressure-volume loop.

### Applicable modes/breath types

All modes and breath types in invasive therapy, except neonate. Not executable when circuit obstruction, circuit disconnect, or high baseline pressure alarm is active. Not executable within 60 seconds of the completion of a Low Flow P/V Maneuver.

### Ranges

Start PEEP: Between 0 and current set PEEP level

Flow rate: 0.5 to 6.0 L/min

$P_{LIMIT}$ : 10 to 60 cmH<sub>2</sub>O

$V_{LIMIT}$ : Adults 100 to 3000 mL, Pediatric 20 to 1000 mL

### To prepare for a Low Flow P/V Maneuver:

- Select the **Low Flow P/V maneuver** app button.
- Set the required control settings within the maneuver app dialogue box.
  - Press **Start PEEP** button and input the desired starting PEEP for the maneuver.
  - Press **Flow rate** button and input desired constant flow for the maneuver.
  - Press  **$P_{LIMIT}$**  button and input desired maximum permissible inspiratory pressure limit.
  - Press  **$V_{LIMIT}$**  button and input desired maximum permissible volume delivery limit.

### To perform a Low Flow P/V Maneuver:

- Press and hold the **HOLD TO START** button to begin the maneuver.
- The ventilator will change PEEP to the *Start PEEP* setting level and deliver a constant flow according to the *Flow Rate* setting, until one of the following conditions is met:
  - The clinician released the **HOLD TO START** button
  - The set  $P_{LIMIT}$  is met
  - The set  $V_{LIMIT}$  is met
  - The set *High Paw* alarm is violated
- When the maneuver is completed, the ventilator will resume normal ventilation.

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